

The Dornier **SideFocus Fiber**

- *High power fiber design optimized for prostate vaporization*
- *Unique side-fire fiber for prostate applications*
- *Supported by different scope manufacturers*

- 
- *High visibility aiming beam*
 - *Indicator line provides accuracy for precise targeting*
 - *Fiber tip positioning control*

Headquarters

Dornier MedTech GmbH
Argelsrieder Feld 7
D-82234 Wessling
Germany

America

Dornier MedTech America, Inc.
1155 Roberts Boulevard
Suite 100
Kennesaw, GA 30144
USA
Phone: +1-770 426-1315
Fax: +1-770 426 6115
Email: info@dornier.com

Asia/Pacific

Dornier MedTech Asia Pte. Ltd.
2 Corporation Road
#05-11/12 Corporation Place
Singapore 618494
Phone: +65-6665-5788
Fax: +65-6665 6711
Email: infoasia@dornier.com

Japan

Dornier MedTech Japan
Meguro Estate Building
3-8-5 Kami-Ohsaki
Shinagawa-ku, Tokyo 141-0021
Japan
Phone: +81-3-3280 3550
Fax: +81-3-3280 3555
Email: info@domedtech.co.jp

Europe

Dornier MedTech Europe GmbH
Argelsrieder Feld 7
D-82234 Wessling
Germany
Phone: +49-8153-888 625
Fax: +49-8153-888 444
Email: infoeurope@dornier.com

Dornier MedTech Italia s.r.l.
Via Arrigo Cavaglieri N. 26
00173 Rome
Italy
Phone: +39-06-72 35 289
Fax: +39-06-72 35 302

Dornier MedTech Espana
Pedro Y Pons 9-11
Edificio Masters
E-08034 Barcelona
Spain
Phone: +34-93-203 9316
Fax: +34-93-205 3814

Dornier MedTech France SARL
156 rue des Castors
74330 Epagny
France
Phone: +33-450-22 18 94
Fax: +33-450-24 24 63

Dornier MedTech Europe GmbH
Moscow Representative Office
Uliza Mytnaya 3
Office 21
119049 Moskau
Russia
Phone: +7-495-230 29 94
Fax: +7-495-739 51 33



CE 1275



The Dornier Medilas D UroBeam

High Power Therapy for BPH

 **LPS.**



 **Dornier MedTech**

WE'RE ALL ABOUT PEOPLE.

The Dornier **Medilas D UroBeam**

The Dornier **Medilas D UroBeam** is a high powered laser optimized for the ablation of prostate tissue in patients with benign prostatic hyperplasia (BPH). This minimally invasive procedure is performed in a hospital or surgery center under spinal or general anesthesia. The advanced technology of the **UroBeam**, combined with 250 watt power and absorption characteristics of the 940 nm

wavelength, excels in making this laser ideal for reducing the symptoms of BPH with minimal side effects.

The **UroBeam**'s state of the art technology and ergonomic design make this an exceptional value. The color touch screen and operator programmable keys provide ease of use, while the small footprint of the laser allows for easy transport.




Patient and Physician Benefits

- *Reduced hospital stay*
- *Minimal catheter time*
- *Optimal hemostasis, minimal bleeding*
- *Treatment of patients with anti-coagulant medication*


Technology Advantages

- *Optimal wavelength to treat tissue independent of hemoglobin content*
- *Green aiming beam for enhanced visibility*
- *Easy to use color touch panel*
- *Small footprint*


Results



1
Inner Coagulation Zone




2
Outer Coagulation Zone



3
Non-Affected Tissue




Canine Specimen: 2cm long from top to bottom



Before

Patient Treatment



After

Clinical pictures: Courtesy of Dr. Thomas Bayer, Kempten, Germany
Microscopic images: Courtesy of Drs. Michael Seitz, Munich, and Thomas Bayer, Kempten, Germany.

Fibertom™ Mode

Utilizing an opto-electronic control system, the proprietary Fibertom mode, in conjunction with Dornier's bare fiber, ensures precise, consistent cutting performance by controlling the fiber tip temperature within defined limits. This feature provides a well controlled application of energy while at the same time minimizing the carbonization on the fiber tip, eliminating tissue attachment to the fiber, and preventing damage to the fiber tip.

Lightguide Protection System

Patient safety is increased with Dornier's patented Lightguide Protection System (LPS). This feature detects light emissions originating from excess heat at the bare fiber tip preventing fiber destruction.



Treatment options for the Medilas D UroBeam

- *Prostate Tissue Ablation for BPH Treatment*
- *Condylomata Acuminata*
- *Kidney Tumors*
- *Penile Carcinoma*
- *Bladder Tumors*
- *Stenoses of Urethra or Ureter*