


AUTOLOGOUS CELLULAR REJUVENATION

INNOVATIVE PROGRAM
FOR SKIN CELL
BIO-STIMULATION



RegenACR® treatment:

Revitalizes and regenerates dermis




Smooths out wrinkles and folds

Induces new collagen production



Regen Lab SA
En Budron B2
CH-1052 Le Mont-sur-Lausanne
SWITZERLAND

T. +41 21 864 38 00
F. + 41 21 864 38 01
info@regenlab.com
www.regenlab.com





Product	Reference	Description	Set/Component
<p>RegenACR® Classic 1 Safety-Lok™ Butterfly needle, 1 Collection holder, 1 RegenBCT tube*, 1x 2 ml ampoule of 10% calcium gluconate, 2 Luer-Lok™ syringes, 1 ml (1) and 5 ml (10), 1 red transfer needle, 3 x 30 G needles for injection, 1 x 4 mm, 27 G needle for mesotherapy, 1 x 80 mm transfer cannula, 1 Transfer device.</p>	<p>R-ACR C1/B</p>	<p>This minimal invasive skin treatment is designed to withdraw blood from the patient and prepare 4 ml of Autologous Platelet Rich Plasma (A-PRP). It's recommended for mesotherapy, fine wrinkles and alopecia treatments.</p>	
<p>RegenACR® Extra 1 Safety-Lok™ Butterfly needle , 1 Collection holder, 2 RegenBCT tubes*, 1 x 2 ml ampoule of 10% calcium gluconate, 2 Luer-Lok™ syringes, 1 ml (1) and 5 ml (10), 1 Red transfer needle, 3 x 30 G needles for injection, 1 x 4 mm, 27 G needle for mesotherapy, 1 x 80 mm transfer cannula, 1 Transfer device, 1 Compressed paper mask.</p>	<p>R-ACR C2/B</p>	<p>RegenACR® Extra is designed to prepare 2 x 4 ml of Autologous Platelet Rich Plasma (A-PRP). It's recommended for expanded mesotherapy as well as fine wrinkles & alopecia treatments. Also includes a compressed paper mask for facial application once poured with A-PRP.</p>	
<p>RegenACR® Plus 1 Safety-Lok™ Butterfly needle, 1 Collection holder, 2 RegenBCT tubes*, 1 RegenATS tube* , 1 Transfer device, 1 x 1 ml Luer-Lok™ syringe, 1 x 80 mm transfer cannula, 3 x 30 G needles for injection, 1 x 4 mm, 27 G needle for mesotherapy, 1 Red transfer needle, 2 x 5 ml Luer-Lok™ syringes.</p>	<p>R-ACR C/BA</p>	<p>RegenACR®Plus is designed to prepare 2 x 4 ml Autologous Platelet Rich Plasma (A-PRP) and 1 x 4 ml autologous thrombin serum (ATS). It's recommended for mesotherapy, fine wrinkle and fold augmentation treatments.</p>	
<p>*Tubes capacity = 8ml</p>			

Product	Reference	Description	Set/Component
<p>RegenACR® Cellular Mask Advanced</p> <p>1 Safety-Lok™ butterfly needle, 1 collection holder, 1 RegenBCT tube*, 1 Cytolnat® Centella cream, 1 spatula.</p>	<p>R-ACR CMC/B</p>	<p>This non-invasive skin treatment includes a special cream which can be applied on face, décolleté and hands after being mixed with plasma. Recommended post-laser, post peeling and post-mesotherapy.</p>	
<p>*Tubes capacity = 8ml</p>			

RegenACR

APPLICATIONS

point of care, in **10** minutes



Product	Reference	Description	Set/Component
RegenCentrilab	RCL-80-2B	<p>Tabletop centrifuge, max speed: 4000 r/min. Rotor imbalance protective device. Door switch & mechanical safety lock.</p> 	
RegenCentrilab	RCL-80-2C	<p>Tabletop centrifuge, max speed: 4000 r/min. Rotor imbalance protective device. Door switch & mechanical safety lock.</p> 	
PRP-Centri	RGL-PRP-C	<p>Tabletop centrifuge, max speed: 4000 r/min. Rotor imbalance protective device. Door switch & mechanical safety lock.</p> 	